

HPX

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW: 51kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human HPX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HX;

Catalog #: A5603

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 3263

Isotype: IgG

Swiss Prot: P02790

Purity: Affinity purification

For research use only.

Background:

This gene encodes a plasma glycoprotein that binds heme with high affinity. The encoded protein is an acute phase protein that transports heme from the plasma to the liver and may be involved in protecting cells from oxidative stress.

To place an order, please [Click HERE](#).