

ECE1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:87kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ECE1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ECE;

Background:

The protein encoded by this gene is involved in proteolytic processing of endothelin precursors to biologically active peptides. Mutations in this gene are associated with Hirschsprung disease, cardiac defects and autonomic dysfunction. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.

To place an order, please [Click HERE](#).

Catalog #:A5638

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1889

Isotype:IgG

Swiss Prot:P42892

Purity:Affinity purification

For research use only.