

F2R

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:47kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human F2R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TR; HTR; CF2R; PAR1; PAR-1;

Background:

Coagulation factor II receptor is a 7-transmembrane receptor involved in the regulation of thrombotic response. Proteolytic cleavage leads to the activation of the receptor. F2R is a G-protein coupled receptor family member.

To place an order, please [Click HERE](#).

Catalog #:A5641

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2149

Isotype:IgG

Swiss Prot:P25116

Purity:Affinity purification

For research use only.