

TNFRSF13B

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:32kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human TNFRSF13B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CVID; RYZN; TACI; CD267; CVID2; IGAD2; TNFRSF14B;

Background:

The protein encoded by this gene is a lymphocyte-specific member of the tumor necrosis factor (TNF) receptor superfamily. It interacts with calcium-modulator and cyclophilin ligand (CAML). The protein induces activation of the transcription factors NFAT, AP1, and NF-kappa-B and plays a crucial role in humoral immunity by interacting with a TNF ligand. This gene is located within the Smith-Magenis syndrome region on chromosome 17.

To place an order, please [Click HERE](#).

Catalog #:A5658

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:23495

Isotype:IgG

Swiss Prot:O14836

Purity:Affinity purification

For research use only.