

DOK1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:52kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human DOK1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P62DOK;

Catalog #:A5687

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1796

Isotype:IgG

Swiss Prot:Q99704

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is part of a signal transduction pathway downstream of receptor tyrosine kinases. The encoded protein is a scaffold protein that helps form a platform for the assembly of multiprotein signaling complexes. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).