

## FOXP3

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**47kDa

**Observed MW:**Refer to figures

**Immunogen:**

Recombinant protein of human FOXP3

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

k

**Synonym:**

JM2; AIID; IPEX; PIDX; XPID; DIETER;

**Catalog #:**A5706

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**50943

**Isotype:**IgG

**Swiss Prot:**Q9BZS1

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a member of the forkhead/winged-helix family of transcriptional regulators. Defects in this gene are the cause of immunodeficiency polyendocrinopathy, enteropathy, X-linked syndrome (IPEX), also known as X-linked autoimmunity-immunodeficiency syndrome. Alternatively spliced transcript variants encoding different isoforms have been identified.

*To place an order, please [Click HERE](#).*