

## CDX1

---

**Reactivity:** Human Mouse

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:1000

**Calculated MW:** 28kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human CDX1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

k

**Synonym:**

Homeobox protein CDX-1; Caudal-type homeobox protein 1;

**Catalog #:** A5712

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 1044

**Isotype:** IgG

**Swiss Prot:** P47902

**Purity:** Affinity purification

For research use only.

**Background:**

This gene is a member of the caudal-related homeobox transcription factor gene family. The encoded DNA-binding protein regulates intestine-specific gene expression and enterocyte differentiation. It has been shown to induce expression of the intestinal alkaline phosphatase gene, and inhibit beta-catenin/T-cell factor transcriptional activity.

*To place an order, please [Click HERE](#).*