

RUVBL1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:20 - 1:50

Calculated MW: 50kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RUVBL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RVB1; TIH1; ECP54; TIP49; INO80H; NMP238; PONTIN; TIP49A; Pontin52;

To place an order, please [Click HERE](#).

Catalog #: A5723

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8607

Isotype: IgG

Swiss Prot: Q9Y265

Purity: Affinity purification

For research use only.