

CDK5

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:33kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CDK5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bi

Synonym:

PSSALRE;

To place an order, please [Click HERE](#).

Catalog #:A5730

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1020

Isotype:IgG

Swiss Prot:Q00535

Purity:Affinity purification

For research use only.