

## PRPF31

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**55kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human PRPF31

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

RP11; PRP31; SNRNP61; NY-BR-99;

**Catalog #:**A5732

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**26121

**Isotype:**IgG

**Swiss Prot:**Q8WWY3

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes a component of the spliceosome complex and is one of several retinitis pigmentosa-causing genes. When the gene product is added to the spliceosome complex, activation occurs.

*To place an order, please [Click HERE](#).*