

GABRR1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:55kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human GABRR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

g

Catalog #:A5745

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2569

Isotype:IgG

Swiss Prot:P24046

Purity:Affinity purification

For research use only.

Background:

GABA is the major inhibitory neurotransmitter in the mammalian brain where it acts at GABA receptors, which are ligand-gated chloride channels. GABRR1 is a member of the rho subunit family. Several transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).