

PCGF6

Reactivity:Human Mouse

Tested applications:WB IF

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:100

Calculated MW:39kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human PCGF6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MBLR; RNF134;

Background:

The protein encoded by this gene contains a RING finger motif, which is most closely related to those of polycomb group (PcG) proteins RNF110/MEL-18 and BMI1. PcG proteins are known to form protein complexes and function as transcription repressors. This protein has been shown to interact with some PcG proteins and act as a transcription repressor. The activity of this protein is found to be regulated by cell cycle dependent phosphorylation. Alternatively spliced transcript variants encoding different isoforms have been identified.

To place an order, please [Click HERE](#).

Catalog #:A5760

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:84108

Isotype:IgG

Swiss Prot:Q9BYE7

Purity:Affinity purification

For research use only.