

NR5A2

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:61kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NR5A2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

k

Synonym:

B1F; CPF; FTF; B1F2; LRH1; LRH-1; FTZ-F1; hB1F-2; FTZ-F1beta;

To place an order, please [Click HERE](#).

Catalog #:A5766

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:2494

Isotype:IgG

Swiss Prot:O00482

Purity:Affinity purification

For research use only.