

## CDKAL1

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

**Calculated MW:** 65kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human CDKAL1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Catalog #:** A5771

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 54901

**Isotype:** IgG

**Swiss Prot:** Q5VV42

**Purity:** Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.

**To place an order, please [Click HERE](#).**