

CDKAL1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:65kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CDKAL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A5771

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:54901

Isotype:IgG

Swiss Prot:Q5VV42

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes.

To place an order, please [Click HERE](#).