

## GSTO2

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200

**Calculated MW:**28kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human GSTO2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

GSTO 2-2; bA127L20.1;

**Catalog #:**A5774

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**119391

**Isotype:**IgG

**Swiss Prot:**Q9H4Y5

**Purity:**Affinity purification

For research use only.

**Background:**

The protein encoded by this gene is an omega class glutathione S-transferase (GST). GSTs are involved in the metabolism of xenobiotics and carcinogens. Four transcript variants encoding different isoforms have been found for this gene.

*To place an order, please [Click HERE](#).*