

TOPBP1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:100 - 1:200 IF 1:50 - 1:200

Calculated MW:170kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human TOPBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

TOP2BP1;

Catalog #:A5781

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11073

Isotype:IgG

Swiss Prot:Q92547

Purity:Affinity purification

For research use only.

Background:

This gene encodes a binding protein which interacts with the C-terminal region of topoisomerase II beta. This interaction suggests a supportive role for this protein in the catalytic reactions of topoisomerase II beta through transient breakages of DNA strands.

To place an order, please [Click HERE](#).