

CD22

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:10 - 1:100

Calculated MW:95kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CD22

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SIGLEC2; SIGLEC-2;

To place an order, please [Click HERE](#).

Catalog #:A5790

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:933

Isotype:IgG

Swiss Prot:P20273

Purity:Affinity purification

For research use only.