

SERPINA6

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:45kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human SERPINA6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CBG;

Catalog #:A5803

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:866

Isotype:IgG

Swiss Prot:P08185

Purity:Affinity purification

For research use only.

Background:

This gene encodes an alpha-globulin protein with corticosteroid-binding properties. This is the major transport protein for glucorticoids and progestins in the blood of most vertebrates. The gene localizes to a chromosomal region containing several closely related serine protease inhibitors which may have evolved by duplication events.

To place an order, please [Click HERE](#).