

CXXC1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF ChIP

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 ChIP 1:20 - 1:100

Calculated MW: 75kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CXXC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

fk

Synonym:

CFP1; CGBP; SPP1; PCCX1; PHF18; hCGBP; ZCGPC1; HsT2645; 2410002I16Rik;
5830420C16Rik;

Catalog #: A5814

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 30827

Isotype: IgG

Swiss Prot: Q9P0U4

Purity: Affinity purification

For research use only.

Background:

Proteins that contain a CXXC motif within their DNA-binding domain, such as CXXC1, recognize CpG sequences and regulate gene expression (Carlone and Skalnik, 2001)

To place an order, please [Click HERE](#).