

RMDN3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 52kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant Protein of human RMDN3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

b

Synonym:

RMD3; RMD-3; FAM82C; FAM82A2; ptpip51;

Catalog #: A5820

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55177

Isotype: IgG

Swiss Prot: Q96TC7

Purity: Affinity purification

For research use only.

Background:

RMDN3 (regulator of microtubule dynamics 3) is a protein-coding gene. Diseases associated with RMDN3 include cerebritis, and noma. GO annotations related to this gene include protein binding. An important paralog of this gene is RMDN1.

To place an order, please [Click HERE](#).