

FIS1

Reactivity:Human Mouse Rat

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW:16kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human FIS1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bh

Synonym:

TTC11; CGI-135;

Catalog #:A5821

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:51024

Isotype:IgG

Swiss Prot:Q9Y3D6

Purity:Affinity purification

For research use only.

Background:

The balance between fission and fusion regulates the morphology of mitochondria. TTC11 is a component of a mitochondrial complex that promotes mitochondrial fission (James et al., 2003

To place an order, please [Click HERE](#).