

## KDM4B

---

**Reactivity:**Human Mouse Rat

**Tested applications:**WB IHC IF

**Recommended Dilution:**WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

**Calculated MW:**121kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human KDM4B

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Synonym:**

JMJD2B; TDRD14B;

*To place an order, please [Click HERE](#).*

**Catalog #:**A5824

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**23030

**Isotype:**IgG

**Swiss Prot:**O94953

**Purity:**Affinity purification

For research use only.