

## HELLS

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IHC IF IP CHIP

**Recommended Dilution:** WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200 IP 1:50 - 1:200

ChIP 1:20 - 1:100

**Calculated MW:** 97kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human HELLS

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

f

**Synonym:**

LSH; PASG; SMARCA6; Nbla10143;

**Catalog #:** A5831

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 3070

**Isotype:** IgG

**Swiss Prot:** Q9NRZ9

**Purity:** Affinity purification

For research use only.

**Background:**

This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.

**To place an order, please [Click HERE](#).**