

RHOT1

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 71kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human RHOT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bh

Synonym:

ARHT1; MIRO1; MIRO-1;

To place an order, please [Click HERE](#).

Catalog #: A5838

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 55288

Isotype: IgG

Swiss Prot: Q8IXI2

Purity: Affinity purification

For research use only.