

RTN4R

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:50kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human RTN4R

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NGR; NOGOR;

Background:

This gene encodes the receptor for reticulon 4, oligodendrocyte myelin glycoprotein and myelin-associated glycoprotein. This receptor mediates axonal growth inhibition and may play a role in regulating axonal regeneration and plasticity in the adult central nervous system.

To place an order, please [Click HERE](#).

Catalog #:A5847

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:65078

Isotype:IgG

Swiss Prot:Q9BZR6

Purity:Affinity purification

For research use only.