

TFRC

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:84kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human TFRC

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

T9; TR; TFR; p90; CD71; TFR1; TRFR;

To place an order, please [Click HERE](#).

Catalog #:A5865

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:7037

Isotype:IgG

Swiss Prot:P02786

Purity:Affinity purification

For research use only.