

Pan DiMethyl-lysine

Reactivity:Human Mouse Rat Other (Wide Range)

Tested applications:WB IF IP

Recommended Dilution:WB 1:500 - 1:2000 IF 1:50 - 1:100 IP 1:50 - 1:100

Observed MW:Refer to Figures

Immunogen:

A synthetic methylated peptide

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Pan methyllysine;

To place an order, please [Click HERE](#).

Catalog #:A5870

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Isotype:IgG

Purity:Affinity purification

For research use only.