

HNRNPR

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:100 - 1:200

Calculated MW: 70kDa

Observed MW: Refer to Figures

Immunogen:

A synthetic peptide of human HNRNPR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HNRPR; hnRNP-R;

Catalog #: A5883

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10236

Isotype: IgG

Swiss Prot: O43390

Purity: Affinity purification

For research use only.

Background:

This gene encodes an RNA-binding protein that is a member of the spliceosome C complex, which functions in pre-mRNA processing and transport. The encoded protein also promotes transcription at the c-fos gene. Alternative splicing results in multiple transcript variants. There are pseudogenes for this gene on chromosomes 4, 11, and 10.

To place an order, please [Click HERE](#).