

## CHD2

**Reactivity:** Human Mouse Rat

**Tested applications:** WB IF

**Recommended Dilution:** WB 1:500 - 1:2000 IF 1:20 - 1:50

**Calculated MW:** 211kDa

**Observed MW:** Refer to Figures

**Immunogen:**

Recombinant protein of human CHD2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

fj

**Synonym:**

EEOC;

**Catalog #:** A5895

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 1106

**Isotype:** IgG

**Swiss Prot:** O14647

**Purity:** Affinity purification

For research use only.

**Background:**

The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.

**To place an order, please [Click HERE](#).**