

NOLC1

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:1000 - 1:2000 IHC 1:50 - 1:200

Calculated MW:73kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human NOLC1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

P130; NOPP130; NOPP140; NS5ATP13;

To place an order, please [Click HERE](#).

Catalog #:A5899

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9221

Isotype:IgG

Swiss Prot:Q14978

Purity:Affinity purification

For research use only.