

EIF1AX

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:200 - 1:2000 IHC 1:100 - 1:200

Calculated MW:16kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human EIF1AX

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

EIF1A; EIF4C; eIF-1A; eIF-4C; EIF1AP1;

Background:

This gene encodes an essential eukaryotic translation initiation factor. The protein is required for the binding of the 43S complex (a 40S subunit, eIF2/GTP/Met-tRNAi and eIF3) to the 5' end of capped RNA.

To place an order, please [Click HERE](#).

Catalog #:A5917

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1964

Isotype:IgG

Swiss Prot:P47813

Purity:Affinity purification

For research use only.