

## PWP2

---

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:200 - 1:2000

**Calculated MW:**102kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant protein of human PWP2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

UTP1; PWP2H; EHOC-17;

*To place an order, please [Click HERE](#).*

**Catalog #:**A5929

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**5822

**Isotype:**IgG

**Swiss Prot:**Q15269

**Purity:**Affinity purification

For research use only.