

SRPR

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 70kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human SRPR

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

DP; Sralpha;

Catalog #: A5939

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6734

Isotype: IgG

Swiss Prot: P08240

Purity: Affinity purification

For research use only.

Background:

The gene encodes a subunit of the endoplasmic reticulum signal recognition particle receptor that, in conjunction with the signal recognition particle, is involved in the targeting and translocation of signal sequence tagged secretory and membrane proteins across the endoplasmic reticulum. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).