

FXR1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:2000

Calculated MW:70kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human FXR1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bj

Synonym:

FXR1P;

Catalog #:A5942

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8087

Isotype:IgG

Swiss Prot:P51114

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is an RNA binding protein that interacts with the functionally-similar proteins FMR1 and FXR2. These proteins shuttle between the nucleus and cytoplasm and associate with polyribosomes, predominantly with the 60S ribosomal subunit.

Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).