

RBM6

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:129kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human RBM6

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

3G2; g16; DEF3; DEF-3; HLC-11; NY-LU-12;

To place an order, please [Click HERE](#).

Catalog #:A5954

Antibody Type:

Monoclonal Antibody

Species:Mouse

Gene ID:10180

Isotype:IgG

Swiss Prot:P78332

Purity:Affinity purification

For research use only.