

STRAP

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:1000 - 1:2000

Calculated MW:38kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human STRAP

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

MAWD; PT-WD; UNRIP;

To place an order, please [Click HERE](#).

Catalog #:A5964

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:11171

Isotype:IgG

Swiss Prot:Q9Y3F4

Purity:Affinity purification

For research use only.