

## ERAL1

---

**Reactivity:**Human

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:2000

**Calculated MW:**48kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human ERAL1

**Storage Buffer:**

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

ERA; H-ERA; ERAL1A; HERA-A; HERA-B;

**To place an order, please [Click HERE](#).**

**Catalog #:**A5974

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**26284

**Isotype:**IgG

**Swiss Prot:**O75616

**Purity:**Affinity purification

For research use only.