

GRWD1

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:49kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human GRWD1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

CDW4; GRWD; RRB1; WDR28;

Catalog #:A6002

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:83743

Isotype:IgG

Swiss Prot:Q9BQ67

Purity:Affinity purification

For research use only.

Background:

This gene encodes a glutamate-rich protein that contains five WD-repeat motifs. The encoded protein may play a critical role in ribosome biogenesis and may also play a role in histone methylation through interactions with CUL4-DDB1 ubiquitin E3 ligase.

To place an order, please [Click HERE](#).