

ALYREF

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:26kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human ALYREF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

ALY; BEF; REF; THOC4; ALY/REF;

Catalog #:A6010

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10189

Isotype:IgG

Swiss Prot:Q86V81

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene is a heat stable, nuclear protein and functions as a molecular chaperone. It is thought to regulate dimerization, DNA binding, and transcriptional activity of basic region-leucine zipper (bZIP) proteins.

To place an order, please [Click HERE](#).