

## AQR

---

**Reactivity:**Human Mouse

**Tested applications:**WB ICC IP

**Recommended Dilution:**WB 1:500 - 1:1000 ICC 1:50 - 1:200 IP 1:20 - 1:50

**Calculated MW:**171kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human AQR

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

IBP160; fSAP164;

*To place an order, please [Click HERE](#).*

**Catalog #:**A6011

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**9716

**Isotype:**IgG

**Swiss Prot:**O60306

**Purity:**Affinity purification

For research use only.