

CELF2

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:500

Calculated MW: 54kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CELF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ETR3; ETR-3; NAPOR; CUGBP2; BRUNOL3;

Catalog #: A6016

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10659

Isotype: IgG

Swiss Prot: O95319

Purity: Affinity purification

For research use only.

Background:

Members of the CELF/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing, and translation. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).