

CELF2

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:500

Calculated MW:54kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human CELF2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ETR3; ETR-3; NAPOR; CUGBP2; BRUNOL3;

Background:

Members of the CELF/BRUNOL protein family contain two N-terminal RNA recognition motif (RRM) domains, one C-terminal RRM domain, and a divergent segment of 160-230 aa between the second and third RRM domains. Members of this protein family regulate pre-mRNA alternative splicing and may also be involved in mRNA editing, and translation. Alternative splicing results in multiple transcript variants encoding different isoforms.

To place an order, please [Click HERE](#).

Catalog #:A6016

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10659

Isotype:IgG

Swiss Prot:O95319

Purity:Affinity purification

For research use only.