

G3BP2

Reactivity:Human

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:500

Calculated MW:54kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human G3BP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

To place an order, please [Click HERE](#).

Catalog #:A6026

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:9908

Isotype:IgG

Swiss Prot:Q9UN86

Purity:Affinity purification

For research use only.