

LIN28A

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:23kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human LIN28A

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CSDD1; LIN28; LIN-28; ZCCHC1; lin-28A;

To place an order, please [Click HERE](#).

Catalog #:A6034

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:79727

Isotype:IgG

Swiss Prot:Q9H9Z2

Purity:Affinity purification

For research use only.