

NOVA2

Reactivity:Human Rat

Tested applications:WB IP

Recommended Dilution:WB 1:500 - 1:2000 IP 1:20 - 1:50

Calculated MW:49kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human NOVA2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ANOVA; NOVA3;

To place an order, please [Click HERE](#).

Catalog #:A6039

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:4858

Isotype:IgG

Swiss Prot:Q9UNW9

Purity:Affinity purification

For research use only.