

## PABPN1

**Reactivity:**Human Mouse

**Tested applications:**WB

**Recommended Dilution:**WB 1:500 - 1:1000

**Calculated MW:**32kDa

**Observed MW:**Refer to Figures

**Immunogen:**

A synthetic peptide of human PABPN1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Concentration:**

j

**Synonym:**

OPMD; PAB2; PABII; PABP2; PABP-2;

**Catalog #:**A6041

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**8106

**Isotype:**IgG

**Swiss Prot:**Q86U42

**Purity:**Affinity purification

For research use only.

**Background:**

This gene encodes an abundant nuclear protein that binds with high affinity to nascent poly(A) tails. The protein is required for progressive and efficient polymerization of poly(A) tails at the 3' ends of eukaryotic transcripts and controls the size of the poly(A) tail to about 250 nt. At steady-state, this protein is localized in the nucleus whereas a different poly(A) binding protein is localized in the cytoplasm. This gene contains a GCG trinucleotide repeat at the 5' end of the coding region, and expansion of this repeat from the normal 6 copies to 8-13 copies leads to autosomal dominant oculopharyngeal muscular dystrophy (OPMD) disease. Related pseudogenes have been identified on chromosomes 19 and X. Read-through transcription also exists between this gene and the neighboring upstream BCL2-like 2 (BCL2L2) gene.

*To place an order, please [Click HERE](#).*