

PAIP1

Reactivity:Human Mouse Rat

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:53kDa

Observed MW:Refer to Figures

Immunogen:

Recombinant protein of human PAIP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #:A6042

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10605

Isotype:IgG

Swiss Prot:Q9H074

Purity:Affinity purification

For research use only.

Background:

The protein encoded by this gene interacts with poly(A)-binding protein and with the cap-binding complex eIF4A. It is involved in translational initiation and protein biosynthesis. Overexpression of this gene in COS7 cells stimulates translation. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified.

To place an order, please [Click HERE](#).