

## PIWIL2

---

**Reactivity:**Mouse Rat

**Tested applications:**WB

**Recommended Dilution:**WB 1:1000 - 1:3000

**Calculated MW:**109kDa

**Observed MW:**Refer to Figures

**Immunogen:**

Recombinant Protein of human PIWIL2

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

CT80; HILL; mili; PIWIL1L;

**Background:**

PIWIL2 belongs to the Argonaute family of proteins, which function in development and maintenance of germline stem cells (Sasaki et al., 2003)

*To place an order, please [Click HERE](#).*

**Catalog #:**A6044

**Antibody Type:**

Polyclonal Antibody

**Species:**Rabbit

**Gene ID:**55124

**Isotype:**IgG

**Swiss Prot:**Q8TC59

**Purity:**Affinity purification

For research use only.