

SNRNP70

Reactivity:Human Mouse

Tested applications:WB

Recommended Dilution:WB 1:500 - 1:1000

Calculated MW:51kDa

Observed MW:Refer to Figures

Immunogen:

A synthetic peptide of human SNRNP70

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

RPU1; Snp1; U1AP; U170K; U1RNP; RNPU1Z; SNRP70; U1-70K;

To place an order, please [Click HERE](#).

Catalog #:A6065

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:6625

Isotype:IgG

Swiss Prot:P08621

Purity:Affinity purification

For research use only.